

QHCC

Generated by Doxygen 1.12.0

1 Namespace Index	1
1.1 Package List	1
2 File Index	3
2.1 File List	3
3 Namespace Documentation	5
3.1 demo_QHCC_XGBoost Namespace Reference	5
3.1.1 Detailed Description	5
3.1.2 Function Documentation	5
3.1.2.1 auto_feature_select()	5
3.1.2.2 experiment_redcap_lab()	6
3.1.2.3 experiment_redcap_lab_radiomics()	6
3.1.2.4 main()	7
3.2 extract_QHCC_radiomics_features Namespace Reference	7
3.2.1 Detailed Description	8
3.2.2 Variable Documentation	8
3.2.2.1 extractor	8
3.2.2.2 feature_row	8
3.2.2.3 features	8
3.2.2.4 flat_features	8
3.2.2.5 image	8
3.2.2.6 index	8
3.2.2.7 mask	9
3.2.2.8 mask_dir	9
3.2.2.9 mask_files	9
3.2.2.10 mask_path	9
3.2.2.11 mri_dir	9
3.2.2.12 mri_files	9
3.2.2.13 mri_path	9
3.2.2.14 output_path	9
3.2.2.15 radiomics_data	9
3.2.2.16 temp_mri_path	9
4 File Documentation	11
4.1 demo_QHCC_XGBoost.py File Reference	11
4.2 extract_QHCC_radiomics_features.py File Reference	11

Chapter 1

Namespace Index

1.1 Package List

Here are the packages with brief descriptions (if available):

demo_QHCC_XGBoost	
The demo file for doing classification on QHCC datasets	5
extract_QHCC_radiomics_features	
Extracting radiomics features from medical images with given masks	7

Chapter 2

File Index

2.1 File List

Here is a list of all files with brief descriptions:

demo_QHCC_XGBoost.py	11
extract_QHCC_radiomics_features.py	11

Chapter 3

Namespace Documentation

3.1 demo_QHCC_XGBoost Namespace Reference

The demo file for doing classification on QHCC datasets.

Functions

- [auto_feature_select](#) (X0, y0, k=1)
Feature selection to select the top k features with the highest mutual information values Input:
- [experiment_redcap_lab](#) (df, seed)
Experiments on prediction with tabular data Input:
- [experiment_redcap_lab_radiomics](#) (df, df_1, df_2, ids_train=None, ids_test=None, seed=100)
Experiments on prediction with both tabular data and radiomics features Input:
- [main](#) ()
The main function for performing experiments and printing out cross validation results.

3.1.1 Detailed Description

The demo file for doing classification on QHCC datasets.

Version

v1.0

3.1.2 Function Documentation

3.1.2.1 auto_feature_select()

```
demo_QHCC_XGBoost.auto_feature_select (  
    X0,  
    y0,  
    k = 1)
```

Feature selection to select the top k features with the highest mutual information values Input:

Parameters

<i>X0</i>	dataframe, the tabular data
<i>y0</i>	dataframe, the label data
<i>k</i>	integer, the number of feature elements to be selected Output:

Returns

`selected_features`: the indexes of the selected features

3.1.2.2 experiment_redcap_lab()

```
demo_QHCC_XGBoost.experiment_redcap_lab (
    df,
    seed)
```

Experiments on prediction with tabular data Input:

Parameters

<i>df</i>	dataframe, the tabular data
<i>seed</i>	random seed Output:

Returns

`model`: the trained model
`ids_train`: the patient IDs for the training data
`ids_test`: the patient IDs for the test data
`acc`: accuracy
`auc`: AUC (area under the curve) score

3.1.2.3 experiment_redcap_lab_radiomics()

```
demo_QHCC_XGBoost.experiment_redcap_lab_radiomics (
    df,
    df_1,
    df_2,
    ids_train = None,
    ids_test = None,
    seed = 100)
```

Experiments on prediction with both tabular data and radiomics features Input:

Parameters

<i>df</i>	dataframe, the tabular data
<i>df_1</i>	dataframe, the radiomics features extracted from MRI images
<i>df_2</i>	dataframe, the radiomics features extracted from CT images
<i>ids_train</i>	the patient IDs for the training data
<i>ids_test</i>	the patient IDs for the test data
<i>seed</i>	random seed Output:

Returns

model: the trained model
 ids_train: the patient IDs for the training data
 ids_test: the patient IDs for the test data
 acc: accuracy
 auc: AUC (area under the curve) score

3.1.2.4 main()

```
demo_QHCC_XGBoost.main ()
```

The main function for performing experiments and printing out cross validation results.

3.2 extract_QHCC_radiomics_features Namespace Reference

Extracting radiomics features from medical images with given masks.

Variables

- `extractor` = `featureextractor.RadiomicsFeatureExtractor()`
Initialize the PyRadiomics feature extractor.
- `str mri_dir` = `"/data/L1_files"`
File paths: MRI image and Mask file directories.
- `str mask_dir` = `"/data/L1_files_mask"`
- `mri_files` = `sorted([f for f in os.listdir(mri_dir) if f.endswith(".nii")])`
Get lists of MRI and Mask files default file format: ".nii".
- `mask_files` = `sorted([f for f in os.listdir(mask_dir) if f.endswith(".nii")])`
- `radiomics_data` = `pd.DataFrame()`
Initialize a DataFrame to store the features.
- `mri_path` = `os.path.join(mri_dir, mri_file)`
Batch processing.
- `mask_path` = `os.path.join(mask_dir, mask_file)`
- `image` = `sitk.ReadImage(mri_path)`
Load MRI image and Mask.
- `mask` = `sitk.ReadImage(mask_path)`
- `temp_mri_path` = `os.path.join(mri_dir, f"temp_{mri_file}")`
Save the adjusted MRI image.
- `features` = `extractor.execute(temp_mri_path, mask_path)`
Extract features.
- `dict flat_features` = `{key: value for key, value in features.items()}`
Flatten features and convert them to a DataFrame row
Use the filename as the ID.
- `feature_row` = `pd.DataFrame([flat_features], index=[mri_file.split('.')[0]])`
- `str output_path` = `"radiomics_features.csv"`
Save Radiomics features to CSV.
- `index`

3.2.1 Detailed Description

Extracting radiomics features from medical images with given masks.

Version

v1.0

3.2.2 Variable Documentation

3.2.2.1 extractor

```
extract_QHCC_radiomics_features.extractor = featureextractor.RadiomicsFeatureExtractor()
```

Initialize the PyRadiomics feature extractor.

3.2.2.2 feature_row

```
extract_QHCC_radiomics_features.feature_row = pd.DataFrame([flat_features], index=[mri_file.↵  
split('.') [0]])
```

3.2.2.3 features

```
extract_QHCC_radiomics_features.features = extractor.execute(temp_mri_path, mask_path)
```

Extract features.

3.2.2.4 flat_features

```
dict extract_QHCC_radiomics_features.flat_features = {key: value for key, value in features.↵  
items() }
```

Flatten features and convert them to a DataFrame row
Use the filename as the ID.

3.2.2.5 image

```
extract_QHCC_radiomics_features.image = sitk.ReadImage(mri_path)
```

Load MRI image and Mask.

3.2.2.6 index

```
extract_QHCC_radiomics_features.index
```

3.2.2.7 mask

```
extract_QHCC_radiomics_features.mask = sitk.ReadImage(mask_path)
```

3.2.2.8 mask_dir

```
str extract_QHCC_radiomics_features.mask_dir = "./data/L1_files_mask"
```

3.2.2.9 mask_files

```
extract_QHCC_radiomics_features.mask_files = sorted([f for f in os.listdir(mask_dir) if f.↵↵endswith(".nii")])
```

3.2.2.10 mask_path

```
extract_QHCC_radiomics_features.mask_path = os.path.join(mask_dir, mask_file)
```

3.2.2.11 mri_dir

```
str extract_QHCC_radiomics_features.mri_dir = "./data/L1_files"
```

File paths: MRI image and Mask file directories.

3.2.2.12 mri_files

```
extract_QHCC_radiomics_features.mri_files = sorted([f for f in os.listdir(mri_dir) if f.↵↵endswith(".nii")])
```

Get lists of MRI and Mask files default file format: ".nii".

3.2.2.13 mri_path

```
extract_QHCC_radiomics_features.mri_path = os.path.join(mri_dir, mri_file)
```

Batch processing.

3.2.2.14 output_path

```
extract_QHCC_radiomics_features.output_path = "radiomics_features.csv"
```

Save Radiomics features to CSV.

3.2.2.15 radiomics_data

```
extract_QHCC_radiomics_features.radiomics_data = pd.DataFrame()
```

Initialize a DataFrame to store the features.

Add to the main DataFrame.

3.2.2.16 temp_mri_path

```
extract_QHCC_radiomics_features.temp_mri_path = os.path.join(mri_dir, f"temp_{mri_file}")
```

Save the adjusted MRI image.

Chapter 4

File Documentation

4.1 demo_QHCC_XGBoost.py File Reference

Namespaces

- namespace [demo_QHCC_XGBoost](#)
The demo file for doing classification on QHCC datasets.

Functions

- [demo_QHCC_XGBoost.auto_feature_select](#) (X0, y0, k=1)
Feature selection to select the top k features with the highest mutual information values Input:
- [demo_QHCC_XGBoost.experiment_redcap_lab](#) (df, seed)
Experiments on prediction with tabular data Input:
- [demo_QHCC_XGBoost.experiment_redcap_lab_radiomics](#) (df, df_1, df_2, ids_train=None, ids_test=None, seed=100)
Experiments on prediction with both tabular data and radiomics features Input:
- [demo_QHCC_XGBoost.main](#) ()
The main function for performing experiments and printing out cross validation results.

4.2 extract_QHCC_radiomics_features.py File Reference

Namespaces

- namespace [extract_QHCC_radiomics_features](#)
Extracting radiomics features from medical images with given masks.

Variables

- `extract_QHCC_radiomics_features.extractor` = `featureextractor.RadiomicsFeatureExtractor()`
Initialize the PyRadiomics feature extractor.
- `str extract_QHCC_radiomics_features.mri_dir` = `"/data/L1_files"`
File paths: MRI image and Mask file directories.
- `str extract_QHCC_radiomics_features.mask_dir` = `"/data/L1_files_mask"`
- `extract_QHCC_radiomics_features.mri_files` = `sorted([f for f in os.listdir(mri_dir) if f.endswith(".nii")])`
Get lists of MRI and Mask files default file format: ".nii".
- `extract_QHCC_radiomics_features.mask_files` = `sorted([f for f in os.listdir(mask_dir) if f.endswith(".nii")])`
- `extract_QHCC_radiomics_features.radiomics_data` = `pd.DataFrame()`
Initialize a DataFrame to store the features.
- `extract_QHCC_radiomics_features.mri_path` = `os.path.join(mri_dir, mri_file)`
Batch processing.
- `extract_QHCC_radiomics_features.mask_path` = `os.path.join(mask_dir, mask_file)`
- `extract_QHCC_radiomics_features.image` = `sitk.ReadImage(mri_path)`
Load MRI image and Mask.
- `extract_QHCC_radiomics_features.mask` = `sitk.ReadImage(mask_path)`
- `extract_QHCC_radiomics_features.temp_mri_path` = `os.path.join(mri_dir, f"temp_{mri_file}")`
Save the adjusted MRI image.
- `extract_QHCC_radiomics_features.features` = `extractor.execute(temp_mri_path, mask_path)`
Extract features.
- `dict extract_QHCC_radiomics_features.flat_features` = `{key: value for key, value in features.items()}`
Flatten features and convert them to a DataFrame row
Use the filename as the ID.
- `extract_QHCC_radiomics_features.feature_row` = `pd.DataFrame([flat_features], index=[mri_file.split('.')[0]])`
- `str extract_QHCC_radiomics_features.output_path` = `"radiomics_features.csv"`
Save Radiomics features to CSV.
- `extract_QHCC_radiomics_features.index`